

Omni - CTR

Capsular Tension Ring



Specification



Capsular Tension Ring		
Material : CLEAR PMMA / BLUE PMMA		
Model	Size	Thickness
CTR 100	10mm	0.18 (±0.02) mm
CTR 110	11mm	0.18 (±0.02) mm
CTR 120	12mm	0.18 (±0.02) mm
CTR 130	13mm	0.18 (±0.02) mm
CTR 140	14mm	0.18 (±0.02) mm

Description:

- One Piece PMMA ring
- High flexibility due to special design
- ETO sterilized
- MATERIAL - PMMA CQUV

Features:

- Circular expansion and stabilization of the periphery of the capsular bag
- Reduced risk of capsular fibrosis
- Reduced risk of capsular shrinking
- Safe IOL centration in eyes with zonular dehiscence

Omni

CTR INJECTOR

The Omni CTR can be comfortably inserted in the capsular bag in a very controlled manner during the surgery with the help of Omni CTR injector system, which is extremely easy to use.



Rated **NSIC-CRISIL SE 1A**

An ISO 9001 & 13485
Certified Company

CE
0434

REGD. OFFICE: _____
OMNI LENS PVT.LTD. 5, "SAMRUDDHI", OPP. SAKAR-III, NAVRANGPURA, AHMEDABAD-380 014. (India)
Ph : +91-79-27544031-32-34 **Fax :** +91-79-27544035 **E-mail :** info@omnilens.in www.omnilens.in

